

Required Forms

1. Cardiology Test Requisition Form

2. Specimen Requirements and Shipping/Handling:

Standard Operating Procedure for blood collection:

Volume: Plastic Whole Blood tube with spray-coated K₂EDTA containing 4 mL of whole sterile blood

Labeling: label with the patient reference, date of collection and hospital reference.

Shipping: performed at room temperature in a padded envelope. Ships immediately for receive within 48h of blood withdrawn (Monday-Thursday).

Company: Any express service (DHL, FEDEX, ...)

Rejection: Blood specimens will be rejected if frozen, hemolyzed or clotted.

Ship to:

HeartGenetics, Genetics and Biotechnology S.A.

A/C Prof. Alexandra Fernandes or Prof. Susana Santos

Taguspark, Parque de Ciência e Tecnologia, Edifício Inovação II, 421, 2740-122 Oeiras

Portugal